

YS6060 Benchtop Grating Spectrophotometer

YS6060 is a benchtop grating spectrophotometer which is developed by our independently with proprietary intellectual property. YS6060 has many features, like 7 inches TFT capacitive touch screen display, full illuminants, reflective d/8 and transmissive d/0 geometry(including or excluding UV). With very stable and precise color measurement, large storage and powerful PC software all makes YS6060 ideal for color analysis within R&D and laboratory environments.



7 inches TFT Capacitive
Touch-screen



USB/Bluetooth



Wavelength range
360nm — 780nm



Built-in camera
locating



PRODUCT FEATURES


- 1.High configuration of hardware: 7 inches TFT Color Capacitive Touch–screen Display; Concave Grating.
- 2.Double Array 256 Image Element CMOS Sensor; Long life–span stable LED, UV LED and xenon lamp.
- 3.With reflective and transmissive spectrum, accurate Lab value, good to calculate color formula and do precise color transmission.
- 4.Auto identify measuring aperture. Freely switchable between 4 measuring apertures: $\phi 25.4\text{mm}/8\text{mm}/4\text{mm}$. Users also can customize apertures.
- 5.Built–in temperature sensor to monitor and compensate the measuring temperature to ensure the measurement more precision.
- 6.Wavelength range 360nm – 780nm. Built–in 400nm cut off/460nm cut off (only xenon lamp edition), more professional in UV measurement.
- 7.Independent light source detector, continuously monitor the condition of light sources to ensure the light source reliable.
- 8.Multiple measurement modes: Quality Management Mode, Sample Mode; Meet more users’ requirement.
- 9.Multiple accessories, sample holders, fixation clamp, suitable to more working condition.
- 10.Big capacity data storage, for 20000 pieces test result.
- 11.Built–in camera locating.
- 12.More powerful extended functions at the PC software.

PRODUCT HIGHLIGHTS



Plastic Sample

Auto–identify the apertures, 25.4/8/4mm three apertures freely switch. Special aperture can be customised.



Transmissive Sample

Transmittance Measurement: Adopt D/0° Geometry, conform to ISO, CIE, ASTM and DIN standard.



Liquid, Powder, Solid Sample Measurement



Measuring Apertures

Auto–identify the apertures, 25.4/15/8/4mm three apertures freely switch. Special aperture can be customised.

APPLICATION INDUSTRY

YS6060 benchtop spectrophotometer is designed for accurate analysis and transmission of laboratory color and appearance of solids, liquids, pastes, powders, pills, and granules in a more comprehensive, streamlined process. YS6060 is widely used in paints, inks, textiles, garments, printing and dyeing, printing etc industries for color transfer and quality control.



Automobile



Leather



Plastics




Paint



Food stuff



Laboratory




Others

YS6060Specification

Optical Geometry:Reflectance: d/8 (SCI&SCE; Include UV/Exclude UV) Transmittance: d/0 (SCI&SCE; Include UV/Exclude UV)		
Comply to standard: CIE No.15, GB/T 3978, GB 2893, GB/T 18833, ISO7724/1, ASTM E1164, DIN5033 Teil7		
Integrating Sphere Size:Φ154mm		
Light Source Device:360nm-780nm Combined LED Lamp, 400nm cut-off, 420nm cut-off, UV Lamp		
Sensor:256 Image Element Double Array CMOS Image Sensor		
Light Source Device:360nm-780nm Combined LED Lamp, 400nm cut-off, 420nm cut-off, UV Lamp		
Wavelength Pitch:10nm		
Semiband Width: 5nm		
Reflectance Range:0~200%		
Measuring Aperture:Reflective: Φ30mm/Φ25.4mm, Φ18mm/Φ15mm, Φ10mm/Φ8mm, Φ6mm/Φ4mm;Transmissive: Φ30mm/Φ25.4mm		
Specular Component:Reflectance: SCI&SCE Transmittance: SCI&SCE		
Color Space:CIE LAB,XYZ,Yxy,LCh,CIE LUV,Musell,s-RGB,HunterLab,βxy,DIN		
Color Difference Formula:ΔE ab,ΔE uv,ΔE 94,ΔE cmc(2:1),ΔE cmc(1:1),ΔE 00, DINΔE99,ΔE(Hunter)		
Colorimetric Index:WI (ASTM E313, CIE/ISO, AATCC, Hunter),YI (ASTM D1925, ASTM 313), MI (Metamerism Index), Staining Fastness, Color Fastness, Color Strength, Opacity,Gardner Index, APHA/Hazen/Pt-Co Index, 555 Index , Saybolt		
Observer Angle:2°&10°		
Illuminants:D65A,C,D50,D55,D75,F1,F2,F3,F4,F5,F6,F7,F8,F9,F10,F11,F12, CWF,DLF,TL83,TL84,TPL5,U30		
Displayed Data:Spectrogram/Values, Chromaticity Values, Color Difference Values/Graph, Pass/Fail Result, Color Offset		
Measurement time: About 2.4 s (simultaneous measurement SCI / SCE about 5 s)		
Repeatability	Spectral reflectance:	Φ25.4mm/SCI Standard deviation within 0.04%(400 nm to 700 nm: within 0.04%)
	Chromaticity value:	Φ25.4mm/SCI Standard deviation within ΔE*ab 0.01 (When a white calibration plate is measured 30 times at 5 second intervals after white calibration)
	Spectral Transmittance:	Φ25.4mm/SCI Standard deviation within 0.05% (400 nm to 700 nm: within 0.04%)
	Chromaticity value:	Φ25.4mm/SCI Standard deviation within ΔE*ab 0.02 (When a white calibration plate is measured 30 times at 5 second intervals after white calibration)
Inter-instrument Error:Φ25.4mm/SCI Within ΔE*ab 0.12(Average for 12 BCRA Series II color tiles)		
Size:370X300X200mm		
Weight:9.6kg		
Power Supply:DC 24V, 3A Power Adapter		
Light Source Device Life:5 years, more than 3 million times measurements.		
Screen:7" TFT Capacitive Screen-touch Display		
Data Port:USB & Bluetooth		
Data Storage Capacity:Standard: 5000 Pcs; Sample: 40000 Pcs.(One PCS can include both SCI and SCE)		
Language:English & Chinese		
Standard Accessory:White and Black Calibration Board, Checking Green Board, Sample Holder, Φ4mm, Φ8mm, Φ15mm, Φ25.4mm Aperture,Power Adapter, USB Cable, User Guide, PC Software, Culture dish , Micro Aperture (4mm) Transmission Test Fixture Assembly		
Optional Accessory: Micro-printers, instrument inverted test fixtures, petri dishes		

GUANGDONG THREENH TECHNOLOGY CO., LTD.



★CONTACT US

-  web:www.3nh.com
-  Email:3nh@3nh.com
-  Tel:0086-020-82880288
-  Add: 6-8th floors, Building B33, Low Carbon Headquarters Park, Xincheng Road No.400, Zengcheng District, Guangzhou, Guangdong Province, China